

## Performance comparison of Phusion Plus and Phusion DNA polymerases

The world leader in serving science



#### Advantages of Thermo Scientific<sup>™</sup> Phusion<sup>™</sup> Plus DNA Polymerase over Phusion<sup>™</sup> DNA polymerases



#### Added benefits:

- No more annealing temperature (T<sub>m</sub>) calculation—uses universal annealing temperature for all primers
- Fewer **PCR runs**—can co-cycle targets of different lengths




#### Improved benefits:


- Improved PCR sequence accuracy—enzyme fidelity increases to >100x that of *Taq* enzyme
- Easier detection of low-abundance targets—higher PCR sensitivity
- Better results with GC-rich sequences—new GC enhancer included
- Higher tolerance to **PCR inhibitors**—reaction buffer reformulated








#### thermofisher.com/phusionplus




#### **Added benefits**

## **No more T<sub>m</sub> calculation**

#### **Thermo Fisher** S C I E N T I F I C






**PCR cycling under two annealing conditions.** 12 targets with varying calculated annealing temperatures (indicated above each lane) were amplified from 50 ng of human genomic DNA (gDNA), following a universal annealing temperature of 60°C (left), or the annealing temperatures calculated with the  $\underline{T}_{\underline{m}}$  calculator (right). The molecular weight marker is  $\underline{Thermo Scientific}^{\underline{m}}$ ZipRuler<sup>TM</sup> Express DNA Ladder 2.


Phusion Plus DNA Polymerase provides convenience and flexibility in PCR preparation

#### **Fewer PCR runs**



Using one cycling protocol


Using different cycling protocols



**PCR co-cycling.** Five targets of different lengths were amplified from human gDNA using a universal cycling protocol for all targets (up to 7.5 kb), with the extension time of the longest amplicon (3 min 45 sec for 7.5 kb) (left) or following separate cycling protocols with a different extension time calculated for each target (9 sec for 0.3 kb, 21 sec for 0.7 kb, 60 sec for 2 kb, 2 min for 3.9 kb, 3 min 45 sec for 7.5 kb) (right). The molecular weight marker is ZipRuler Express DNA Ladder 2.


The universal annealing feature of Phusion Plus DNA Polymerase also allows a **universal cycling protocol**, helping you to circumvent multiple PCR runs and save time. One annealing temperature (60°C) and one extension time based on the longest amplicon can be used for targets of different lengths—i.e., **co-cycling different targets** on the same block—without compromising PCR yields and specificity.

Different PCR targets can be co-cycled using Phusion Plus DNA Polymerase



#### **Improved benefits**

### Improved PCR sequence accuracy






**Tech note:** Learn more about enzyme fidelity measurement using NGS with molecular barcodes.

**Fidelity of high-fidelity polymerases measured relative to** *Taq* **DNA polymerase.** Error rates were determined by next-generation sequencing (NGS) using molecular barcodes, then normalized to that of *Taq* DNA polymerase.

Phusion Plus DNA Polymerase's fidelity is >100x that of *Taq* DNA polymerase, giving you more confidence in PCR sequence accuracy


#### **Easier detection of low-abundance targets**



**High sensitivity of Phusion Plus DNA Polymerase.** A 0.5 kb target was amplified from different amounts of human genomic (gDNA). The molecular weight marker is the <u>Thermo Scientific<sup>™</sup> GeneRuler<sup>™</sup> 1 kb Plus DNA Ladder</u>.

Phusion Plus DNA Polymerase can detect targets from as little as 0.08 ng of human gDNA


#### **Better results with GC-rich sequences**



Efficient GC-rich amplification of Phusion Plus DNA Polymerase. 16 targets with high GC content (their percentages indicated) were amplified from 50 ng of human gDNA. The molecular weight marker is <u>Thermo Scientific<sup>™</sup> ZipRuler<sup>™</sup> Express</u> DNA Ladder 2.


Phusion Plus DNA Polymerase includes the new Phusion<sup>™</sup> GC Enhancer in the package for more efficient amplification of sequences with >65% GC content

#### **Higher tolerance to PCR inhibitors**




High inhibitor tolerance of Phusion Plus DNA Polymerase. A 2 kb target was amplified from 50 ng of human gDNA. The reaction mixtures contained at a final concentration of: 1—no inhibitor, 2—humic acid (0.5  $\mu$ g/mL), 3—hemin (2.5  $\mu$ M), or 4—xylan (250  $\mu$ g/mL). The molecular weight marker is the <u>GeneRuler 1 kb Plus DNA Ladder</u>.

Phusion Plus DNA Polymerase better tolerates inhibitors from plants (e.g., xylan), soil (e.g., humic acid), and blood (e.g., hemin) during PCR



#### **Comparable performance**

## High yields and specificity




**High yields and specificity using Phusion Plus DNA Polymerase.** 0.3–7.5 kb DNA targets were amplified from 100 ng of human gDNA. The molecular weight marker is the <u>GeneRuler 1 kb Plus</u> <u>DNA Ladder</u>.

Similar PCR specificity and yields can be achieved with Phusion Plus DNA Polymerase

#### Thermo Fisher

### **Benchtop stability of assembled reactions**




Non-hot-start DNA polymerases are not recommended for lab automation setup because assembled reactions are not stable at room temperature for an extended period. DNA Polymerase

**Phusion High-Fidelity** 

**Benchtop stability of Phusion Plus DNA Polymerase.** A 0.5 kb target was amplified from 50 ng of human gDNA. Assembled PCR reactions were loaded immediately onto a thermal cycler (0 hr) or set at room temperature for 24 hr before cycling (24 hr). The molecular weight marker is the <u>GeneRuler 1 kb Plus DNA Ladder</u>.

PCR reactions assembled with Phusion Plus DNA Polymerase are stable at room temperature for up to 24 hr, enabling lab automation setup



## Summary

### Summary: comparison of Phusion Plus, Phusion High-Fidelity, and Phusion Hot Start II DNA polymerases





**Thermo Fisher** 

	Phusion Plus DNA Polymerase	Phusion High-Fidelity DNA Polymerase	Phusion Hot Start II DNA Polymerase
Fidelity (vs. <i>Taq</i> DNA polymerase)	>100x	50x	50x
Hot-start modification (Affibody molecule-mediated)	Yes	No	Yes
Universal annealing temperature (no $T_m$ calculator needed)	Yes	No	No
Universal cycling protocol (co-cycling targets of different length)	Yes	No	No
PCR sensitivity	+++	++	++
GC-rich amplification	+++ (New GC enhancer)	++	++
Inhibitor tolerance	+++	++	++
PCR yields and specificity	+++	+++	+++
Benchtop stability (of assembled reactions)	Up to 24 hr	N/A	Up to 24 hr
Stand-alone and master mix formats	Yes	Yes	Yes

Find your upgrade at thermofisher.com/phusionplus

# Thank you

For Research Use Only. Not for use in diagnostic procedures. © 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. COL114090 0321