

ARE YOU BEING SERVED?

THE BENEFITS OF DUAL DATA
ACQUISITION IN A ROUTINE DIOXIN
LAB



In cooperation with

ThermoFisher
SCIENTIFIC

DSP systems **MIURA**
The missing link in your laboratory!


NofaLab

- Strategically positioned in the harbour
- In the vicinity of some of the largest plant oil refineries in Europe

Europe





Rotterdam area




START OF DIOXINS AT NOFALAB

- 2008 1st DFS and 1st generation cleanup system (FMS)
- 2010 2nd DFS and two 1st generation cleanup system
- 2013 both DFS systems equipped with Dual Data Acquisition (D.D.A.) (1st generation)
- 2015 3rd DFS installed with 2nd generation D.D.A. and 2 (2nd generation) Cleanup systems



Samples per month



Kind of samples:

1. Vegetable/Animal Oils, fats and Fatty Acids
2. Human food (eg. Milk, eggs, fish, meat, seeds)
3. Animal feed (eg. compound feed, trace materials)
4. Ect.



HOW DO WE DO IT?

Extraction



Concentrate



Measurement

Clean-up



HOW DO WE DO IT



THIS IS HOW WE DID IT: FROM 72 HOURS REPORTING TIME TO 24 AND 8 HOURS

- 2008-2013 without DDA
 - 72 hours standard reporting time
 - 10 samples a day 24 hours reporting time


- 2013-20xx with DDA:
 - 75% of the samples within 24hours

- 2014-20xx with DDA and GO-HT
 - 8 samples within 8 hours Dioxins and PCBs



8 HOUR SERVICE (ONLY FOR OIL TYPE SAMPLES)

- | | |
|---|----------------|
| ○ Weight in sample ($\pm 1,5$ g of fat) | 2 min |
| ○ Spike with ^{13}C labeled std. | 1 min |
| ○ Load sample onto the column | 5 min |
| ○ Attach all columns | 2 min |
| ○ Load column set to the system | 1 min |
| ○ Start program | 73 min |
| ○ Concentrate | 40 min |
| ○ Inject in GC-HRMS (D.D.A.) | 51 min* |
| ○ Evaluating and reporting | 5 min |
| Total: | 180 min |



* After the first injection all other samples take 26 minutes.



STRATEGY FOR 40 SAMPLES A DAY (POSSIBLE)?

- Preparing the first batch/run most time consuming.
- When first batch will be processed by cleanup system, the 2nd batch can already be prepared.
You have 73 minutes time to finish.
- Working time 7am till 5 pm
- 10 hours / 1,52 hours = 6,58 sequences a day.
- Everyday a blanc system and recovery (1 for oil type samples and 1 for solid type samples) is done
- Staff: 3/4 men working on the extraction/clean-up
2 men do the evaluating/reporting of the data



HOW TO MANAGE

- Average weekly number of samples 160
- On top at least 15 QA samples + 15 Process blanks
- Often 5 samples re-analysis (confirmation)
- Total number of samples injected 400 fractions
- On top at least 15 dielectric standards + 15 DL and NDL- PCB standards
- Total number of injections is 430
- 60 min nonpolar column → run-time for each fraction 55 min
- In case of two single GC-HRMS instruments 197,08 hrs required per instrument


Solution

Lack of > 29,8 hrs per week

ion Module



HOW DO WE MANAGE GC-HRMS DUAL ACQUISITION



Dual acquisition = 1,8x improvement of efficiency
 Every 60 minutes 2 samples are measured



TYPICAL SEQUENCE NOFALAB BV

Week26_diox [Open] - Sequence Setup - Home Page

File Edit Change Actions View GoTo Help

Status	Acquisition Queue	File Name	Path	Inst Meth	Position
1		CS1_190624_GC1	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	1
2		CS1_190624_GC2	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	1
3		CS2_190624_GC1	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	2
4		CS2_190624_GC2	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	2
5		CS2_190624_GC1	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	2
6		CS3_190624_GC2	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	3
7		CS3_190624_GC1	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	3
8		CS4_190624_GC2	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	4
9		CS4_190624_GC1	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	4
10		CS5_190624_GC2	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	5
11		Blanco01	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	6
12		Blanco02	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	7
13		blanco03	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	6
14		Blanco04	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	7
15		blanco05	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	6
16		blanco06	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	7
17		2019030630_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	115
18		2019030631_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	116
19		2019030691_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	117
20		2019029974_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	118
21		2019030460_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	119
22		2019030632_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	120
23		2019030633_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	121
24		2019029983_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	122
25		2019030591_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	123
26		2019030696_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	124
27		2019030693_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	125
28		2019030695_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	126
29		2019030692_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	127
30		2019030694_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	128
31		2019030590_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	129
32		2019030564_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	130
33		2019030594_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	131
34		CS1_V_190625_GC2	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	1



Status	Acquisition Queue	File Name	Path	Inst Meth	Position
		34 CS1_V_190625_GC2	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	1
		35 CS1_V_190625_GC1	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	1
		36 RecSer_190624	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	137
		37 Rec190624	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	140
		38 BlancoSys_SER_190624	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	138
		39 BlancoSys_190624	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	139
		40 2019030589_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	132
		41 2019030254_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	139
		42 2019030608_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	139
		43 2019030595_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	139
		44 2019030893_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	141
		45 2019030372_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	142
		46 2019030364_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	143
		47 2019030757_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	144
		48 2019030572_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	136
		49 2019030848_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	145
		50 2019030662_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	146
		51 2019030661_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	147
		52 2019031127_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	148
		53 2019030981_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	149
		54 2019030955_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	150
		55 CS2_V_190625_GC1	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	2
		56 CS2_V_190625_GC2	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	2
		57 Blanco07	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	7
		58 Blanco08	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	7
		59 2019030933_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	8
		60 2019030742_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	9
		61 2019030745_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	10
		62 2019030965_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	11
		63 2019030749_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	12
		64 2019028426_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	13
		65 2019030565_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	14
		66 2019030759_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	15
		67 2019031069_d	C:\Xcalibur\Data\2019\Week26	C:\Xcalibur\methods\DD_Dioxins_Furans	16
		68			



ADDITIONAL BENEFITS OF THE DUAL DATA OPTION

- Source will stay cleaner for a longer time (no solvent and early eluting analytes do enter the source anymore).
- With the DualData XL option the sample throughput of a DFS can be almost doubled”
- Fast column change (no need to break the vacuum in the source).



DISADVANTAGES OF THE DUAL DATA OPTION

- Once a sequence is started in Dual Data, no changes in sequence order can be made
- No back-up system
 - When you have downtime, you don't have a second system



FOR THE COMING FUTURE?

- The possibility to include a sample and/or update the sequence queue?!



URGENT Sample

The screenshot shows a software interface with a toolbar at the top and a table below. The toolbar includes various icons for file operations and a green 'UPDATE' button with a circular arrow icon. The table has four columns: File Name, Path, Inst Meth, and Position. The first three rows of the table are highlighted in blue. An orange arrow points from the 'UPDATE' button to the 'URGENT Sample' text. Another orange arrow points from the 'URGENT Sample' text to the 'Sequence Row #13' entry in the file explorer on the left.

	File Name	Path	Inst Meth	Position
1	CS1_190715_GC1	C:\Xcalibur\Data\2019\Week29	C:\Xcalibur\methods\DD_Dioxins_Furans	1
2	CS1_190715_GC2	C:\Xcalibur\Data\2019\Week29	C:\Xcalibur\methods\DD_Dioxins_Furans	1
3	CS1_V_190715	C:\Xcalibur\Data\2019\Week29	C:\Xcalibur\methods\DD_Dioxins_Furans	150
4	CS2_190715	C:\Xcalibur\Data\2019\Week29	C:\Xcalibur\methods\DD_Dioxins_Furans	2
5	CS3_190715	C:\Xcalibur\Data\2019\Week29	C:\Xcalibur\methods\DD_Dioxins_Furans	3
6	CS4_190715	C:\Xcalibur\Data\2019\Week29	C:\Xcalibur\methods\DD_Dioxins_Furans	4
7	CS5_190715	C:\Xcalibur\Data\2019\Week29	C:\Xcalibur\methods\DD_Dioxins_Furans	5
8	Blanco01	C:\Xcalibur\Data\2019\Week29	C:\Xcalibur\methods\DD_Dioxins_Furans	6
9	Blanco02	C:\Xcalibur\Data\2019\Week29	C:\Xcalibur\methods\DD_Dioxins_Furans	7
10	blanco03	C:\Xcalibur\Data\2019\Week29	C:\Xcalibur\methods\DD_Dioxins_Furans	6
11	Blanco04	C:\Xcalibur\Data\2019\Week29	C:\Xcalibur\methods\DD_Dioxins_Furans	7
12	2019034928_d	C:\Xcalibur\Data\2019\Week28	C:\Xcalibur\methods\DD_Dioxins_Furans	69
13	2019034872_d	C:\Xcalibur\Data\2019\Week28	C:\Xcalibur\methods\DD_Dioxins_Furans	73
14	2019035054_d	C:\Xcalibur\Data\2019\Week28	C:\Xcalibur\methods\DD_Dioxins_Furans	75
15	2019035060_d	C:\Xcalibur\Data\2019\Week28	C:\Xcalibur\methods\DD_Dioxins_Furans	76
16	2019034953_d	C:\Xcalibur\Data\2019\Week28	C:\Xcalibur\methods\DD_Dioxins_Furans	77
17	2019034955_d	C:\Xcalibur\Data\2019\Week28	C:\Xcalibur\methods\DD_Dioxins_Furans	78
18	2019034962_d	C:\Xcalibur\Data\2019\Week28	C:\Xcalibur\methods\DD_Dioxins_Furans	79
19	2019035061_d	C:\Xcalibur\Data\2019\Week28	C:\Xcalibur\methods\DD_Dioxins_Furans	80
20	2019034963_d	C:\Xcalibur\Data\2019\Week28	C:\Xcalibur\methods\DD_Dioxins_Furans	81
21	2019035165_d	C:\Xcalibur\Data\2019\Week28	C:\Xcalibur\methods\DD_Dioxins_Furans	82

CONCLUSION

- To increase the throughput in your Laboratory and you don't have deep pockets;

DualData XL is the Solution

- Less downtime, because source will stay cleaner for a longer time.



In cooperation with

ThermoFisher
SCIENTIFIC

DSP systems **MIURA**
The missing link in your laboratory!

NofaLab
laboratories



Thank you for your
attention!

Any

Questions??

